



gh2 Antibody

Product Code	CSB-PA323220XA01OEI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20332
Immunogen	Recombinant Oncorhynchus mykiss (Rainbow trout) (Salmo gairdneri) gh2 protein
Raised In	Rabbit
Species Reactivity	Oncorhynchus mykiss (Rainbow trout) (Salmo gairdneri)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oncorhynchus mykiss (Rainbow trout) (Salmo gairdneri)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gh2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.