



sspA Antibody

Product Code	CSB-PA330114XA01FWU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31784
Immunogen	Recombinant Histophilus somni (Haemophilus somnus) sspA protein
Raised In	Rabbit
Species Reactivity	Histophilus somni (Haemophilus somnus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Histophilus somni (Haemophilus somnus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	sspA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.