



Protein translation factor SUI1 homolog Antibody

Product Code	CSB-PA319632XA01MEZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14021
Immunogen	Recombinant Methanococcus vannielii Protein translation factor SUI1 homolog protein
Raised In	Rabbit
Species Reactivity	Methanococcus vannielii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Methanococcus vannielii
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.