



STR1 Antibody

Product Code	CSB-PA300722XA01FXI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68175
Immunogen	Recombinant Rauvolfia serpentina (Serpentine wood) (Ophioxylon serpentinum) STR1 protein
Raised In	Rabbit
Species Reactivity	Rauvolfia serpentina (Serpentine wood) (Ophioxylon serpentinum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rauvolfia serpentina (Serpentine wood) (Ophioxylon serpentinum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	STR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.