



SOD Antibody

Product Code	CSB-PA177213XA01SXQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01137
Immunogen	Recombinant Schistosoma mansoni (Blood fluke) SOD protein
Raised In	Rabbit
Species Reactivity	Schistosoma mansoni (Blood fluke)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Schistosoma mansoni (Blood fluke)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SOD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.