



sodA Antibody

Product Code	CSB-PA022398XA01LGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28763
Immunogen	Recombinant Listeria ivanovii sodA protein
Raised In	Rabbit
Species Reactivity	Listeria ivanovii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Listeria ivanovii
Lead Time	Made-to-order (14-16 weeks)
Gene Names	sodA