



Rz1 Antibody

Product Code	CSB-PA895450XA01ECC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9XJJ7
Immunogen	Recombinant Enterobacteria phage 933W (Bacteriophage 933W) Rz1 protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage 933W (Bacteriophage 933W)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage 933W (Bacteriophage 933W)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Rz1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.