



sufE Antibody

Product Code	CSB-PA484264XA01E00
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7M0N4
Immunogen	Recombinant Escherichia coli O8 (strain IAI1) sufE protein
Raised In	Rabbit
Species Reactivity	Escherichia coli O8 (strain IAI1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli O8 (strain IAI1)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	sufE