



scalR Antibody

Product Code	CSB-PA527495XA01FNZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O52691
Immunogen	Recombinant Streptomyces caespitosus scalR protein
Raised In	Rabbit
Species Reactivity	Streptomyces caespitosus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces caespitosus
Lead Time	Made-to-order (14-16 weeks)
Gene Names	scalR