



tonB Antibody

Product Code	CSB-PA332651XA01KBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45610
Immunogen	Recombinant Klebsiella pneumoniae tonB protein
Raised In	Rabbit
Species Reactivity	Klebsiella pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Klebsiella pneumoniae
Lead Time	Made-to-order (14-16 weeks)
Target Names	tonB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.