



tamB Antibody

Product Code	CSB-PA108842XA01DTP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D2TN57
Immunogen	Recombinant Citrobacter rodentium (strain ICC168) (Citrobacter freundii biotype 4280) tamB protein
Raised In	Rabbit
Species Reactivity	Citrobacter rodentium (strain ICC168) (Citrobacter freundii biotype 4280)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Citrobacter rodentium (strain ICC168) (Citrobacter freundii biotype 4280)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	tamB