



tgm2 Antibody

Product Code	CSB-PA023462XA01EQR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52181
Immunogen	Recombinant Pagrus major (Red sea bream) (Chrysophrys major) tgm2 protein
Raised In	Rabbit
Species Reactivity	Pagrus major (Red sea bream) (Chrysophrys major)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pagrus major (Red sea bream) (Chrysophrys major)
Lead Time	Made-to-order (14-16 weeks)
Target Names	tgm2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.