



Sporulated oocyst TA4 antigen Antibody

Product Code	CSB-PA318360XA01EJU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13399
Immunogen	Recombinant Eimeria tenella (Coccidian parasite) Sporulated oocyst TA4 antigen protein
Raised In	Rabbit
Species Reactivity	Eimeria tenella (Coccidian parasite)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Eimeria tenella (Coccidian parasite)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.