



pacB Antibody

Product Code	CSB-PA716856XA01EDI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5XLR0
Immunogen	Recombinant Enterobacteria phage P7 (Bacteriophage P7) pacB protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage P7 (Bacteriophage P7)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage P7 (Bacteriophage P7)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pacB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.