



ITEVIIIIR Antibody

Product Code	CSB-PA656420XA01ELJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q38419
Immunogen	Recombinant Enterobacteria phage RB3 (Bacteriophage RB3) ITEVIIIIR protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage RB3 (Bacteriophage RB3)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage RB3 (Bacteriophage RB3)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ITEVIIIIR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.