



Thermophilic serine proteinase Antibody

Product Code	CSB-PA674928XA01FXN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q45670
Immunogen	Recombinant Bacillus sp. (strain AK1) Thermophilic serine proteinase protein
Raised In	Rabbit
Species Reactivity	Bacillus sp. (strain AK1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus sp. (strain AK1)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.