



serS Antibody

Product Code	CSB-PA309082XA01FLF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P95689
Immunogen	Recombinant Staphylococcus aureus (strain NCTC 8325) serS protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus (strain NCTC 8325)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Staphylococcus aureus (strain NCTC 8325)
Lead Time	Made-to-order (14-16 weeks)
Target Names	serS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.