



CYP73A4 Antibody

Product Code	CSB-PA342564XA01DQN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48522
Immunogen	Recombinant Catharanthus roseus (Madagascar periwinkle) (Vinca rosea) CYP73A4 protein
Raised In	Rabbit
Species Reactivity	Catharanthus roseus (Madagascar periwinkle) (Vinca rosea)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Catharanthus roseus (Madagascar periwinkle) (Vinca rosea)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CYP73A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.