



# tetB Antibody

Product Code	CSB-PA320390XA01SPJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14551
Immunogen	Recombinant Streptomyces rimosus tetB protein
Raised In	Rabbit
Species Reactivity	Streptomyces rimosus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces rimosus
Lead Time	Made-to-order (14-16 weeks)
Gene Names	tetB