



tgt Antibody

Product Code	CSB-PA019146XA01EOP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7MPG7
Immunogen	Recombinant Escherichia coli O81 (strain ED1a) tgt protein
Raised In	Rabbit
Species Reactivity	Escherichia coli O81 (strain ED1a)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli O81 (strain ED1a)
Lead Time	Made-to-order (14-16 weeks)
Target Names	tgt
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.