



trxA Antibody

Product Code	CSB-PA025365XA01AHD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80579
Immunogen	Recombinant Alicyclobacillus acidocaldarius (Bacillus acidocaldarius) trxA protein
Raised In	Rabbit
Species Reactivity	Alicyclobacillus acidocaldarius (Bacillus acidocaldarius)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Alicyclobacillus acidocaldarius (Bacillus acidocaldarius)
Lead Time	Made-to-order (14-16 weeks)
Target Names	trxA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.