



sacIR Antibody

Product Code	CSB-PA516234XA01FNS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O31074
Immunogen	Recombinant Streptomyces achromogenes sacIR protein
Raised In	Rabbit
Species Reactivity	Streptomyces achromogenes
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces achromogenes
Lead Time	Made-to-order (14-16 weeks)
Target Names	sacIR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.