



VINST1 Antibody

Product Code	CSB-PA328792XA01VFQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28343
Immunogen	Recombinant Vitis vinifera (Grape) VINST1 protein
Raised In	Rabbit
Species Reactivity	Vitis vinifera (Grape)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vitis vinifera (Grape)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	VINST1