



TAT Antibody

Product Code	CSB-PA878699XA01FTB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9M6F0
Immunogen	Recombinant Taxus cuspidata (Japanese yew) TAT protein
Raised In	Rabbit
Species Reactivity	Taxus cuspidata (Japanese yew)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Taxus cuspidata (Japanese yew)
Lead Time	Made-to-order (14-16 weeks)
Target Names	TAT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.