



gene Antibody

Product Code	CSB-PA149956XA01FRR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O64331
Immunogen	Recombinant Escherichia phage N15 (Bacteriophage N15) gene protein
Raised In	Rabbit
Species Reactivity	Escherichia phage N15 (Bacteriophage N15)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage N15 (Bacteriophage N15)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gene
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.