



Triosephosphate isomerase Antibody

Product Code	CSB-PA024102XA01GFW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36186
Immunogen	Recombinant Giardia intestinalis (Giardia lamblia) Triosephosphate isomerase protein
Raised In	Rabbit
Species Reactivity	Giardia intestinalis (Giardia lamblia)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Giardia intestinalis (Giardia lamblia)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.