



Triosephosphate isomerase, cytosolic Antibody

Product Code	CSB-PA325158XA01CXF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21820
Immunogen	Recombinant Coptis japonica (Japanese goldthread) Triosephosphate isomerase, cytosolic protein
Raised In	Rabbit
Species Reactivity	Coptis japonica (Japanese goldthread)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Coptis japonica (Japanese goldthread)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.