



Triosephosphate isomerase, glycosomal Antibody

Product Code	CSB-PA361369XA01TQV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04789
Immunogen	Recombinant Trypanosoma brucei brucei Triosephosphate isomerase, glycosomal protein
Raised In	Rabbit
Species Reactivity	Trypanosoma brucei brucei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Trypanosoma brucei brucei
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.