



Trehalose phosphorylase Antibody

Product Code	CSB-PA524919XA01FXQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75003
Immunogen	Recombinant Grifola frondosa (Maitake) (Polyporus frondosus) Trehalose phosphorylase protein
Raised In	Rabbit
Species Reactivity	Grifola frondosa (Maitake) (Polyporus frondosus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Grifola frondosa (Maitake) (Polyporus frondosus)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.