



38 Antibody

Product Code	CSB-PA220514XA01FTP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XQC4
Immunogen	Recombinant Escherichia phage T1 (Bacteriophage T1) 38 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage T1 (Bacteriophage T1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage T1 (Bacteriophage T1)
Lead Time	Made-to-order (14-16 weeks)
Target Names	38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.