



Tropomyosin Antibody

Product Code	CSB-PA889038XA01BZJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NAS5
Immunogen	Recombinant Anisakis simplex (Herring worm) Tropomyosin protein
Raised In	Rabbit
Species Reactivity	Anisakis simplex (Herring worm)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Anisakis simplex (Herring worm)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.