



tus Antibody

Product Code	CSB-PA525152XA01KBI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O52715
Immunogen	Recombinant Klebsiella pneumoniae subsp. ozaenae tus protein
Raised In	Rabbit
Species Reactivity	Klebsiella pneumoniae subsp. ozaenae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Klebsiella pneumoniae subsp. ozaenae
Lead Time	Made-to-order (14-16 weeks)
Target Names	tus
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.