



ulaG Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA809927XA01FQR |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8FAJ4 |
| Immunogen | A synthesized peptide derived from Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) ulaG protein |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) |
| Tested Applications | ELISA, Dot blot (guaranteed specific to the immunizing peptide) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) |
| Lead Time | Made-to-order (14-16 weeks) |
| Target Names | ulaG |
| Product Modify | ①. 100ug * KLH-conjugated peptide paired with 100ug * BSA-conjugated peptide (positive control); ②. 1ml*pre-immune serum (used as negative control) for free; ③. Rabbit polyclonal antibodies purified by Antigen Affinity; |
| Product Precision | 1. CUSABIO guarantees the antibody purity > 95% confirmed by SDS-PAGE. 2. CUSABIO guarantees the antibody titer > 1: 64,000 confirmed by ELISA. 3. CUSABIO guarantees to offer Peptide-specific Dot blot. 4. CUSABIO guarantees 100% risk-free: No qualified antibody, no charge. |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |