



TRXH Antibody

Product Code	CSB-PA304452XA01DSJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80028
Immunogen	Recombinant Chlamydomonas reinhardtii (Chlamydomonas smithii) TRXH protein
Raised In	Rabbit
Species Reactivity	Chlamydomonas reinhardtii (Chlamydomonas smithii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Chlamydomonas reinhardtii (Chlamydomonas smithii)
Lead Time	Made-to-order (14-16 weeks)
Target Names	TRXH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.