



# ttcA Antibody

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA486751XA01EOR                                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                 |
| <b>Uniprot No.</b>         | B7LRU9                                                                                        |
| <b>Immunogen</b>           | Recombinant Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) ttcA protein |
| <b>Raised In</b>           | Rabbit                                                                                        |
| <b>Species Reactivity</b>  | Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)                          |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)                                                  |
| <b>Form</b>                | Liquid                                                                                        |
| <b>Conjugate</b>           | Non-conjugated                                                                                |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4              |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                     |
| <b>Isotype</b>             | IgG                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                    |
| <b>Product Type</b>        | Polyclonal Antibody                                                                           |
| <b>Immunogen Species</b>   | Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)                          |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)                                                                   |
| <b>Target Names</b>        | ttcA                                                                                          |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                   |