



Uricase-2 isozyme 1 Antibody

Product Code	CSB-PA356320XA01GGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04670
Immunogen	Recombinant Glycine max (Soybean) (Glycine hispida) Uricase-2 isozyme 1 protein
Raised In	Rabbit
Species Reactivity	Glycine max (Soybean) (Glycine hispida)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.