



uvrB Antibody

Product Code	CSB-PA318011XA01FXT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10125
Immunogen	Recombinant <i>Micrococcus luteus</i> (strain ATCC 4698 / DSM 20030 / JCM 1464 / NBRC 3333 / NCIMB 9278 / NCTC 2665 / VKM Ac-2230) (<i>Micrococcus lysodeikticus</i>) uvrB protein
Raised In	Rabbit
Species Reactivity	<i>Micrococcus luteus</i> (strain ATCC 4698 / DSM 20030 / JCM 1464 / NBRC 3333 / NCIMB 9278 / NCTC 2665 / VKM Ac-2230) (<i>Micrococcus lysodeikticus</i>)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	<i>Micrococcus luteus</i> (strain ATCC 4698 / DSM 20030 / JCM 1464 / NBRC 3333 / NCIMB 9278 / NCTC 2665 / VKM Ac-2230) (<i>Micrococcus lysodeikticus</i>)
Lead Time	Made-to-order (14-16 weeks)
Target Names	uvrB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.