



vdh Antibody

Product Code	CSB-PA687483XA01SNX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53872
Immunogen	Recombinant Streptomyces cinnamonensis vdh protein
Raised In	Rabbit
Species Reactivity	Streptomyces cinnamonensis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces cinnamonensis
Lead Time	Made-to-order (14-16 weeks)
Gene Names	vdh