



uvrA Antibody

Product Code	CSB-PA321129XA01MTO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13567
Immunogen	Recombinant Micrococcus luteus (Micrococcus lysodeikticus) uvrA protein
Raised In	Rabbit
Species Reactivity	Micrococcus luteus (Micrococcus lysodeikticus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Micrococcus luteus (Micrococcus lysodeikticus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	uvrA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.