



uvrD Antibody

Product Code	CSB-PA929889XA01SXB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05311
Immunogen	Recombinant Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720) uvrD protein
Raised In	Rabbit
Species Reactivity	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Lead Time	Made-to-order (14-16 weeks)
Target Names	uvrD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.