



rIIA Antibody

Product Code	CSB-PA356119XA01EDZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03690
Immunogen	Recombinant Enterobacteria phage T4 (Bacteriophage T4) rIIA protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T4 (Bacteriophage T4)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rIIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.