



32 Antibody

Product Code	CSB-PA022706XA01ELN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O21945
Immunogen	Recombinant Enterobacteria phage RB70 (Bacteriophage RB70) 32 protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage RB70 (Bacteriophage RB70)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage RB70 (Bacteriophage RB70)
Lead Time	Made-to-order (14-16 weeks)
Target Names	32
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.