



vioA Antibody

Product Code	CSB-PA147952XA01SZA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6U1I3
Immunogen	Recombinant Shigella dysenteriae vioA protein
Raised In	Rabbit
Species Reactivity	Shigella dysenteriae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Shigella dysenteriae
Lead Time	Made-to-order (14-16 weeks)
Target Names	vioA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.