



xis Antibody

Product Code	CSB-PA303939XA01FWP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68927
Immunogen	Recombinant Escherichia phage HK022 (Bacteriophage HK022) xis protein
Raised In	Rabbit
Species Reactivity	Escherichia phage HK022 (Bacteriophage HK022)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage HK022 (Bacteriophage HK022)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	xis