



virF Antibody

Product Code	CSB-PA313898XA01YAQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C2V5
Immunogen	Recombinant Yersinia enterocolitica virF protein
Raised In	Rabbit
Species Reactivity	Yersinia enterocolitica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Yersinia enterocolitica
Lead Time	Made-to-order (14-16 weeks)
Target Names	virF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.