



hkaC Antibody

Product Code	CSB-PA860175XA01FXV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9AZ38
Immunogen	Recombinant Salmonella phage HK620 (Bacteriophage HK620) hkaC protein
Raised In	Rabbit
Species Reactivity	Salmonella phage HK620 (Bacteriophage HK620)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Salmonella phage HK620 (Bacteriophage HK620)
Lead Time	Made-to-order (14-16 weeks)
Target Names	hkaC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.