



xyn2 Antibody

Product Code	CSB-PA151400XA01HQM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HGE1
Immunogen	Recombinant Humicola grisea var. thermoidea xyn2 protein
Raised In	Rabbit
Species Reactivity	Humicola grisea var. thermoidea
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Humicola grisea var. thermoidea
Lead Time	Made-to-order (14-16 weeks)
Target Names	xyn2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.