



frd.2 Antibody

Product Code	CSB-PA762342XA01FXA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76VG9
Immunogen	Recombinant Enterobacteria phage LZ3 (Bacteriophage LZ3) frd.2 protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage LZ3 (Bacteriophage LZ3)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage LZ3 (Bacteriophage LZ3)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	frd.2