



Snake venom serine protease salmobin Antibody

Product Code	CSB-PA525550XA01GGQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O73800
Immunogen	Recombinant Gloydius halys (Chinese water mocassin) (Agkistrodon halys) Snake venom serine protease salmobin protein
Raised In	Rabbit
Species Reactivity	Gloydius halys (Chinese water mocassin) (Agkistrodon halys)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gloydius halys (Chinese water mocassin) (Agkistrodon halys)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.