



xylR Antibody

Product Code	CSB-PA328199XA01CAL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40981
Immunogen	Recombinant Caldicellulosiruptor sp. (strain Rt8B.4) xylR protein
Raised In	Rabbit
Species Reactivity	Caldicellulosiruptor sp. (strain Rt8B.4)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caldicellulosiruptor sp. (strain Rt8B.4)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	xylR